













Dijkous – Bas cuisse

Maat Taille	Grootste kuitomtrek La plus grande circ. du mollet		Beenlengte (hiel-bilplooi) Longueur de jambe (talon-pli fessier)		Kleuren Code couleur	Ref. Réf.
	 cm		 cm			
1	small	23 - 30	regular	64-80		46335
2	small	23 - 30	long	80-97		46336
3	medium	30 - 38	regular	64-80		46337
4	medium	30 - 38	long	80-97		46338
5	large	38 - 46	regular	64-80		46388
6	large	38 - 46	long	80-97		46389
7	x-large	46 - 53	regular	64-80		77125
8	x-large	46 - 53	long	80-97		77128
9	xx-large	53 - 66	regular	64-80		77130
10	xx-large	53 - 66	long	80-97		77131

Alle therapeutische drukkousen waarvan de maat niet correct is kunnen complicaties tijdens het gebruik veroorzaken

Tout bas de compression thérapeutique dont la taille n'est pas correcte peut entraîner des complications